

Ultrasol AF



Application: According to the Regulation of biocide products: Disinfection and cleaning of medical inventory and all surfaces. According to the MDD (Medical Device Directive) 93/42 / EC: Disinfection and cleaning of non-invasive medical devices product aldehyde.
Use: disinfection of medical devices and surfaces is performed: scrubbing and cleaning or using cleaning machines.

Concentration and exposure time: The bactericidal effect of including MRSA, fungicidal according to EN 1372 / EN 13624 / EN 13697: 0.25% (2.5 g / l) - 15 min.

Composition: 1.68 g didecildimetilamonijumhlid 1.1 g of N- (3-aminopropyl) -N-dodecylpropane-1,3-diamine 3.33 g quaternary ammonium compound is benzyl-C12-16-alkildimetilhlid Contains: nonionic surfactants